

Foral™ 85-E CG Hydrogenated Rosinate

FORAL™ 85-E CG hydrogenated gum rosin is a cosmetic grade resin derived from the esterification of a highly stabilized gum rosin and glycerol. This light amber, thermoplastic resin has excellent resistance to oxidation and discoloration caused by heat and aging. FORAL 85-E CG contributes adhesion to various cosmetic and depilatory wax formulations.

For further information regarding this product please refer to:

Synthomer Adhesive Technologies

eMail: Adhesive.Technologies@Synthomer.com

- Contributes to adhesion
- Derived from a natural, renewable source
- Excellent oxidative and color stability
- Light color
- Low odor
- Wide solubility and compatibility range

Property	Typical Value	Unit	Method ¹
Description, Base Resin	Hydrogenated Gum Rosin		
Ring and Ball Softening Point	74	°C	ASTM E 28
Color, Gardner	3		ASTM D 6166, 50% solids in toluene
Acid Number	7	mg KOH/g	ASTM D 465
Density at 25°C	1.06	kg/dm ³	

¹ internal method based upon the specified norm

Applications

Personal Care Additives

Packaging

FORAL™ 85-E CG hydrogenated rosin is pastillated and packed in polyethylene bags of 25 kg net.

Storage

Due to the thermoplastic behavior, pastillated and flaked resins may fuse, block or lump. This can be accelerated under any of the following conditions: 1) above ambient temperature 2) prolonged storage 3) pressure, e.g., stacking pallets, or a combination of these conditions. This is particularly applicable for low softening point resin grades. In order to maintain the flake or pastille shape, we therefore recommend storing the material in a temperature-controlled area; be careful with stacking material or applying pressure and preventing prolonged storage. It should be noted that lumping does not have a negative impact on the product specifications. Due to the nature of the product, claims regarding lumping cannot be accepted.

Resins are prone to gradual oxidation, some more so than others. This could result in darkening and/or it could have an adverse effect on the solubility of the resin in organic solvents or on its compatibility with polymers. Accordingly, it is recommended that strict control of inventory be observed at all times, taking care that the oldest material is used first.

FORAL™ 85-E CG hydrogenated rosin material will remain within product specification limits for a period of at least twelve months after shipment from Synthomer production facilities in the Netherlands, provided storage conditions outlined in this data sheet are observed. However, as we can neither anticipate the conditions under which the resin is processed nor the end-use applications for which it is used, we recommend that the material be tested upon receipt.

Comments

Properties reported here are typical of average lots. Synthomer makes no representation that the material in any particular shipment will

Foral™ 85-E CG Hydrogenated Rosinate



Revision: 11.04.2022
Page 2 of 2

conform exactly to the values given.